

NeuroLIGHT TREATMENT

(888) 426-8794
121 Woodland Ave, Suite 100
Reno NV 89523

The durable light treatment device was designed for frequent clinical use by clinical professionals. Molded and sealed surface, multi-treatment areas and easy-in-and-out to fit any patient quickly. The design reinforces the cables - to insure long life - and comes with an industry leading 3-year warranty.

DEVICE Features

Easy in and out for patients

15 second setup time for clinics

3-year warranty

Capacitance Sensing Technology - switches device on when foot/hand is in place.

Sealed against moisture - reducing damage and worry

Proprietary Progressive Multi-pulse Modulation® and Progressive Pulsed Width Modulation®

Constant current drivers for uniformity in output

Cableless design for easy housekeeping easy (power supply excepted)

Over-molded cable strain relief - integrated into circuit board eliminating plug failure

Embedded microprocessor unit auto terminates program for safety

Multiple Treatment areas - Plantar aspect of the foot, medial and lateral aspect of calf

Worldwide power supply

Use while seated or supine

Accommodates various foot widths & lengths

Made in the USA

Bilateral or unilateral use

Integrated strapping to secure spatulas proximal to leg, base unit for supine

FEATURES Per Side

Diode count PER SIDE

Plantar aspect, each foot	IR 72, Red 31
Legs (combined, each leg)	IR 54, Red 30
TOTAL	187 each side

Emitted Radiation • mw/cm2

IR	3297
Red	2074
Total	5371

Wavelengths

IR	842nm
Red	631nm

